mediroyal 

# NRX<sup>®</sup> ErixThree Neuro Shoulder


MediRoyal Nordic AB  
Box 7052 · SE-192 07 Sollentuna · Sweden  
+46 8 506 766 00 · [www.mediroyal.se](http://www.mediroyal.se)






English	11	<b>GB</b>
Svenska	12	<b>SE</b>
Suomi	14	<b>FI</b>
Deutsch	15	<b>DE</b>
Français	17	<b>FR</b>
Italiano	18	<b>IT</b>
Español	20	<b>ES</b>
Português	21	<b>PT</b>
Nederlands	23	<b>NL</b>
Ελληνικά	25	<b>EL</b>
Eesti	26	<b>ET</b>
Lietuvos	28	<b>LT</b>
Latvijas	29	<b>LV</b>
Polski	31	<b>PL</b>
Русский	32	<b>RU</b>
Český	34	<b>CZ</b>
Arabic	36	<b>AR</b>
Persian	38	<b>PE</b>
Chinese Traditional	40	<b>CH</b>
Chinese Mandarin	42	<b>GM</b>
Japanese	44	<b>JP</b>
Malay	46	<b>MY</b>
Română	48	<b>RO</b>
Slovenski	49	<b>SL</b>

## English · Instructions

### NRX® ErixThree Neuro Shoulder

#### 1. General information

Read the instruction carefully. The product should be fitted by a medical professional like an occupational therapist, physical therapist, orthopedic technician or a medical doctor. Patients with diabetes, circulatory problems or sensitive skin should consult their doctor before use. If the patient is allergic, check the material content on the textile label.

#### 2. Indications

The product has been designed to provide support to the shoulder joint after neurological disorders, stroke or mild to moderate anterior shoulder subluxations.

#### 3. Contraindications

Patients who haven't received a correct diagnosis from a doctor should consult one before use of the brace. The brace will not prevent severe dislocations of the shoulder joint.

#### 4. Pre Application

To select the right size, take the circumference measurement around the mid upper arm. Use the measurement chart to determine the most suitable size. The brace is fitted with a shoulder strap (A) to provide more abduction support, two elastic hook straps (B) that can be used for additional reinforcement and facilitation of selected muscle groups as well as an extender part (C) for the chest strap. On the front part of the brace there are two pockets with flexible steel splints. They can be adjusted individually to prevent pressure over the clavicle bone. Remove them from the pocket and bend them by hand. The inside of the brace is fitted with a friction pattern to increase proprioception.

#### 5. Application

A t-shirt should be worn under the brace. Open all the hook closures on the brace. Slide the arm sleeve up over the arm and position the brace so that the highest part rests on the trapezius muscle. Close the hook parts on the arm and make sure that the textile tongue protects the skin from the hook closures. Don't over tighten the closures as it might disturb circulation. Grab the chest strap and extend it to retract the shoulder posteriorly. Run the strap under the arm pit and attach the strap on the front part of the brace. If the strap is too short, apply the extender to the hook part.

#### 6. Additional abduction support

The shoulder strap (A) can be used to provide more abduction support. Attach one of the hook parts to the loop area on the chest strap on the opposite shoulder. Stretch the strap and apply it over the trapezius, running below C7 over the cervical area and over to the other side. For more abduction support, apply the other hook part on the lateral part of the arm sleeve. For more posterior retraction of the shoulder, apply the hook part over the anterior side of the arm sleeve.

#### 7. Muscle facilitation

The elastic hook straps can be used to facilitate muscles like Supraspinatus, Infraspinatus or the Deltoid. For the best effect, attach one end of the elastic hook over

the insertion area of the muscle, then stretch it slightly (15–20%) and attach over the muscle, attaching it over the area where the muscle ends. You can also try to run the strap the other way from muscle ending to insertion for a different effect.

## **8. Removal**

Open the hook closures on the brace and remove it. Close all straps when the brace is not in use.

## **9. Material content**

Check the label for material content.

## **10. Caring instructions and precautions**

We recommend that the brace should be washed once a week. Follow the washing instructions on the product, close all straps and use a laundry bag. Remove the splints from the front pockets.

Patients with sensitive skin, diabetes or poor circulation should be observed. Patients with allergies should check the material content on the label before use.

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Več informacij je na voljo na [www.mediroyal.se](http://www.mediroyal.se)

mediroyal 

## Mediroyal – Functional Swedish Design

Product development is a combination of ideas, material, design and know-how. At Mediroyal we have more than 20 years of experience with developing products that provides both the patient and the therapist with enhanced functionality, very day. Swedish design is unique in many ways. Quality, performance, functionality and design are all features that make a difference. That's why we are proud to label our products with "Functional Swedish Design". Because orthoses is our passion.